

# NMR Sample Manager for TopSpin: Quick Reference

**Documentation:** [nmr-samples.github.io/topspin](https://nmr-samples.github.io/topspin)

**Web app:** [nmr-samples.github.io/online](https://nmr-samples.github.io/online)

Record what's in the magnet alongside your experiment data. All fields are optional – annotate as much or as little as you need. Sample metadata is stored as JSON files in your data directories. No data leaves your computer.

## Getting Started

Type `samples` in TopSpin to open the Sample Manager. The window shows a timeline of your experiments, a sample list, and a searchable catalogue.

## Commands

Use these instead of the standard TopSpin commands. They perform the same physical operations and additionally capture sample metadata.

| Command             | Replaces           | What it does  |
|---------------------|--------------------|---|
| <code>ija</code>    | <code>ij</code>    | Injects sample (lift air off), then opens annotation form     |
| <code>eja</code>    | <code>ej</code>    | Ejects sample (lift air on), closes annotation with timestamp |
| <code>sxa 15</code> | <code>sx 15</code> | SampleJet: exchanges sample then opens annotation form        |
| <code>sxa ej</code> | <code>sx ej</code> | SampleJet: ejects and returns to storage, closes annotation   |

## What You Can Record

Everything is optional. Fill in what's relevant for you.

| Section              | Fields  |
|----------------------|---|
| <b>People</b>        | Users, research groups  |
| <b>Sample</b>        | Label, components (names, concentrations, isotopic labelling)                 |
| <b>Buffer</b>        | pH, solvent, chemical shift reference, buffer components                      |
| <b>NMR Tube</b>      | Diameter, tube type (regular, Shigemi, solid-state rotor, ...), sample volume |
| <b>Lab Reference</b> | Sample ID, lab book entry   |
| <b>Notes</b>         | Free-form text  |

## Key Features

**Retrospective annotation:** Right-click experiments in the timeline to create an annotation at any time.

**Duplicate:** select a sample and click Duplicate to copy it as a template (ideal for titration series).

**Catalogue:** search across all your data directories for samples by user, buffer, components, or any text.

**Timeline:** see everything you've run within a directory. Double-click to open an experiment in TopSpin.

**One active sample:** creating a new sample automatically marks the previous one as ejected.

## Feedback Is Welcome!

Please send any feedback and comments to [c.waudby@ucl.ac.uk](mailto:c.waudby@ucl.ac.uk).

Alternatively, raise an issue at <https://www.github.com/nmr-samples/topspin/>.